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Trocken lagern

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## Instructions for use MEDUTEX for groin

Before using the product, please be sure to read all instructions for use. This will help to avoid damage due to usage errors and ensure optimal use.

We reserve the right of error and technical modifications.

All currently valid product information can be found in the download area at [www.johngmbh.com](http://www.johngmbh.com).

### Intended use / indication

MEDUTEX for groin is a compression bandage that prevents post-operative bleeding after medical procedures (e.g. cardiac catheterisation, invasive diagnostics and therapy) by fixing compressors.

MEDUTEX for groin provides compression of punctured femoral arteries and veins after invasive diagnostics and therapy for left and right heart catheterisation.

### Contraindications

Known allergies to one or more of the materials used.

### CE-conformity

John GmbH, as manufacturer, declares under its sole responsibility that MEDUTEX of John GmbH is a medical device in compliance with the applicable requirements of Regulation (EU) 2017/745 of the European Parliament and of the Council.

### Lifetime

The product has components that are subject to natural wear even during normal use. The wear depends very much on the individual intensity of use and the maintenance. Particularly in the case of intensive use, components can reach their wear limit even before the end of the statutory warranty period. In these cases of premature wear due to use, however, there is not automatically a defect in the product.

### Re-use

MEDUTEX for groin: The product can only be reused by a single user and not by further users. 

MEDUTEX Leiste for groin ECO: The product cannot be reused. 

### Disposal

The product must be disposed of in accordance with the relevant national legal requirements.

### Warranty & Liability

Warranty and liability are in accordance with the general terms and conditions (GTC, para. VI) of John GmbH.

### Reporting serious incidents

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the the responsible EU authority. „Serious incident“ means any incident that directly or indirectly led, might have led or might lead to any of the following:

a) the death of a patient, user or other person; b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health; c) a serious public health threat – an event which could result in imminent risk of death, serious deterioration in a person's state of health, or serious illness, that may require prompt remedial action, and that may cause significant morbidity or mortality in humans, or that is unusual or unexpected fort he given place and time. (Regulation (EU) 2017/745)

Keep the products in a dry place

Protect the products from direct sunlight

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## Safety instructions


- ⚠ Follow the care instructions and the storage notes.
- ⚠ MEDUTEX is applicable up to a girth of 130 cm.
- ⚠ Always check the Velcro straps for functional suitability before putting on the product.
- ⚠ Always use help for putting on the bandage by a doctor or a professional.
- ⚠ In case your product has not been put on correctly, do not try to move the product by tugging. The material might be damaged. Open your product again and put it on once more.
- ⚠ Coordinate the wearing time with a doctor or a specialist.
- ⚠ Do not use MEDUTEX in case of intolerance to any of the materials.
- ⚠ Do not make any changes to the product.

## Use (also see the pictures)

1. Put the MEDUTEX around the lying patient, so that the velcro part is right of the body.
2. When the pelvic belt is positioned correctly it has to be closed at the side.
3. After covering the puncture spot in a sterile way, the compression can be put on.
4. Then stretch and fix the thigh belt over the compression to the pelvic belt.
5. Right side.
6. Left side (handling as right side).

## Instructions for care & Material

### MEDUTEX for groin

⚠ Close Velcro fasteners before washing

**Material:** Cotton, rayon, perlon, elastane Velcro fastener and loop layer

### MEDUTEX for groin ECO

⚠ **FOR ONE-WAY USE ONLY**

**Material:** Velours (Foam RG29 on both sides coated with Velours) Velcro fastener and loop layer